

RT ALIGN ROTATE

Radiotherapy Phantoms

Features

- Laser alignment marks
- Fiducial marks for 5x5, 10x10, 15x15, and 20x20 field sizes
- Embedded BBs for phantom edges, centre, and offset locations
- Accurate placement of alignment marks and BBs
- Can be rotated 360° in 45° increments
- Use with or without rotation jig

Benefits

- Isocentre alignment using central embedded BBs
- EPID-based field size and light field coincidence checks
- EPID scaling and movement using embedded BBs
- Use on EPID surface or couch

